Monoprinting By Lea Cioci CPD CPT

Supplies Needed: 2 Clear Acrylic Sheets – (You can use a report cover that can be purchased at an office supply store, as it is easy to sandwich the glossy cardstock between) Dye Reinkers in various colors Glossy Cardstock Stamps and Permanent Ink Brayer

Instructions:

1. Lay the acrylic report cover flat on your surface and open the pages like you were reading a book.

2. Squirt 2-3 colors of the reinkers randomly on the bottom half of the open report cover.

3. Take the glossy cardstock and lay it glossy side down on top of the reinker colors. Since the report cover is a good size, you can lay 2 pieces of glossy cardstock that is cut into 1/4 sheets.

4. Close the top part of the report cover over the glossy cardstock that is glossy side down on the reinker colors.

5. Take the brayer and brayer over the closed report cover.

6. Open the report cover and carefully remove the glossy cardstock. Viola' beautifully blended patches of colors are created much like a Rorschach Blot Test.



7. Stamp on top of your cardstock once it has dried.

